

Otto Bock®

QUALITY FOR LIFE



Skippi

Experience the World

Exceed your limits, increase your independence, conquer the outdoors: Skippi makes it fun to discover something new every day. Skippi's speed and high degree of maneuverability help it support children's urge to move and the pleasure they derive from making discoveries. Mobility expands your horizons. Children learn by being actively involved in their surroundings. Arm and back rest adjustments are variable. In this way, the chair can be set to facilitate „seeing and doing“. And this goes for traveling too: Skippi is easy to disassemble and fits in any small car. What's more: Skippi grows with the child! Adapting its seat depth and width is child's play.

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt · Phone +49 5527 848-1590 · Fax +49 5527 848-1676 · healthcare@ottobock.de · www.ottobock.com



Skippy – Power Wheelchair for children

Options

- Seat tilt mechanism, electric or with gas compression spring
- electrical Backrest angle adjustment
- Seating shell adapter
- Tray
- Control for attendant
- Bumber bar
- Single-panel footrest with aluminum footplate

Technical Data

- Max. speed6 km/h
- Battery capacity 12V, 31Ah C20
- Distance range approx. 30 km
- Max. load capacity50 kg
- Empty weight68 kg
- Max. obstacle height that can be cleared ...50 mm
- Overall width57 cm
- Overall length without legrests70 cm
- Overall length with legrests85 cm
- Turning radius55 cm
- Seat height43 cm
- Seat width26 – 38 cm
- Seat depth30 – 37 cm
- Seat inclination (optional: continuous)-9°/3°/12° (0°– 25°)
- Armrest height16 – 27 cm
- Back height35/40/42/44/46 cm
- Back width30 cm
- Back angle (optional: continuous)0-25°
- Lower leg length11 – 26 cm
- Frame colorsyellow, blue, red



Height adjustable push handles for attendant



Electrical backrest angle adjustment



Seat adaptable to size



Detachable footrests

Otto Bock®

QUALITY FOR LIFE