

INVACARE®

TDX® SI

Power Wheelchair



Yes, you can.®

The Invacare® TDX® SI Wheelchair offers a streamlined approach to fully featured rehab power wheelchairs. With standard features such as Invacare® as MK6i™ electronics, simplified Invacare® SureStep® suspension, stability assist, 22NF battery tray and a top speed of 5 mph, the TDX SI Wheelchair delivers performance in an elegant design. The TDX SI Wheelchair offers several seating options including an adjustable Rehab Seat, Captain's seating and Invacare® Formula™CG Tilt. **With a focus on excellent standard features, the TDX SI Wheelchair is the only name you need to know.**



Invacare® TDX® SI Wheelchair Specifications

Overall dimensions	Adj. ASBA Rehab Seat	Formula™CG Tilt	Van Seat
Base length (w/o front riggings)	35"	35"	35"
Base length (w/front riggings)	38.5"	38.5"	38.5"
Base width n1	24" (w/o joystick)	24" (w/o joystick)	24" (w/o joystick)
Height	35.5" - 39.5"	36.5" - 40.5"	38" w/o headrest 45" w/headrest
Motor package	2-pole	2-pole	2-pole
Seat dimensions			
Width	12" - 22"	12" - 22"	16" - 22"
Depth	12" - 22"	12" - 22"	16" - 20"
Back height	12" - 24"	16" - 24"	18"
with headrest	-	-	25"
Back angle	85° - 105°	85° - 105°	85° - 115°
Seat-to-floor height			
Low	16.5" @ 0°	17.5" @ 0°	21"
Medium	18.5" @ 0°	18.5" @ 0°	N/A
Tall	20.5" @ 0°	19.5" @ 0°	N/A
Seat-to-floor height at 5°	add 0.5" to each	add 1" to each	N/A
Weight			
With batteries	254 lb.	295 lb.	267 lb.
Without batteries	178 lb.	219 lb.	191 lb.
Battery options	22NF	22NF	22NF
Wheel and caster size			
Front/rear casters	6" x 2"	6" x 2"	6" x 2"
Drive wheels	14" x 3"	14" x 3"	14" x 3"
Maximum speed	5 mph	5 mph	5 mph
Incline capability	9°	9°	9°
Weight capacity	300 lb.	300 lb.	300 lb.
Turning radius	23"	23"	23"
Ground clearance	3"	3"	3"
Standard MK6i electronics	MK660/SPJ™+	MK660ACC/SPJ™+ w/PSS	MK660/SPJ™+
Limited warranties			
Base frame	60 months	60 months	60 months
Seat frame	36 months	36 months	36 months
Electronics	12 months	12 months	12 months
Motor/gearbox	12 months	12 months	12 months

n1 90° front riggings

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010
Invacare Corporation
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed to
Invacare Corporation
unless otherwise noted.
Specifications are
subject to change without
notification.

Form No. 07-073
Rev. 07/10

The TDX SI Wheelchair is available in the Group 3 Medicare code and is coded as follows:

TDXSI	(Group 3 Rehab Seat)	K0848
TDXSIV	(Group 3 Captain's Seat)	K0849
TDXSI-CG	(Group 3 Single Power)	K0856

