



Yes, you can.®

## ➤ Invacare® Bed Accessories

Improve safety, improve access, improve control

Invacare Bed Accessories are designed to improve the safety of the bed, enable easier access in and out of the bed or make the bed easier to operate. Choosing the right accessory for a new or existing Invacare bed can be the difference between being able or unable to carry out activities independently.



### Cotside (Side Rail)

#### Designed with patient safety and comfort in mind

Serves as a preventive feature against the patient's falling out of the bed. Can be fitted or removed without the necessity of using any tools. The full length cotsides lock into place and lower to below the mattress. Cotsides come in left and right sides and are ordered separately.

#### Available for:

- Alegio
- Etude Basic (SB385)
- Etude Classic (SB200)
- SB750
- Octave
- Kiwicare



### Cotside (Side Rail) Cover

#### Upgrade Cotside Safety

These padded cotside covers are designed to protect users from hurting themselves on the metal safety cotsides. The covers are made from thick vinyl, that can be easily wiped clean. The covers can be attached to cotsides by velcro fasteners or domed buttons. The velcro version is designed to cover one side of the cotside and can flip over. The domed button version fully encases the cotside on the inside and outside.

#### Available for:

- Alegio
- Etude Basic (SB385)
- Etude Classic (SB200)
- SB750
- Octave
- Kiwicare



### Support Handle

#### Helps make getting in and out of bed easier

The Support Handle enables the use of both hands to sit up on the side of the bed and then also offers a means to maintain balance while rising to a standing position. Support Handles are easy to attach and remove, no tools required. A variety of sizes and lengths are available.

#### Available for:

- Alegio
- Etude Basic (SB385)
- Etude Classic (SB200)
- SB750
- Octave
- Kiwicare



### Lifting Pole

#### Enabling easier transfer in and out of bed

A lifting pole provides overhead support for those who find it difficult to rise out of bed. Designed for raising a person to a sitting position and not lifting the entire weight of the person off the bed. The lifting handle can be easily adjusted in both height and depth.

#### Available for:

- Alegio
- Etude Basic (SB385)
- Etude Classic (SB200)
- SB750
- Octave
- Kiwicare



### Bumper Wheels

#### Protect the bed and the walls

Bumper wheels are a recommend option to protect both the bed and the walls from those little, and sometimes big, knocks and bumps.

#### Available for:

- Alegio
- SB750
- Octave



### Soft Controller

#### Ensures a higher degree of independence and safety

The Invacare Soft Controller is ideal for someone who has limited hand function and/or a degree of visual impairment. Patients with reduced mobility in hands and fingers, will clearly find it easier to operate the Invacare Soft Controller as the buttons are generally bigger and require far less strength to activate as well as to keep down, enabling effortless adjustment of the bed. The coloured buttons lead to an easier understanding of the functionality of the bed.

#### Available for:

- Etude Basic (SB385) – new models
- Etude Classic (SB200) – new models
- SB750 – standard with bed



### Infared Controller

#### No wires to get in the way

Control a Kiwicare using an Infared Controller to activate the functions of the bed. No need for annoying cables that can get caught up, reducing the potential for crushed cords.

#### Available for:

- Kiwicare



### Companion Bed

#### When you need to share your Kiwicare

The Companion bed is a flexi-slat bed that has been designed to attach to a Kiwicare bed. This allows for one side of the bed to operate independently as a medical bed while the other side (the Companion) remains static, giving the overall appearance of a domestic bed.

#### Available for:

- Kiwicare (selected sizes)



### Bed Extensions

#### Lengthen your Etude Bed

Mattress support extension enables the lengthening of the bed to suit the longer person. Beds can be extended by 20cms. Extending the bed is a quick and easy procedure, simply split the bed base in half and insert the extension.

#### Available for:


- Etude Basic (SB385)
- Etude Classic (SB200)

## Important information for ordering Invacare Bed Accessories

- In most cases Accessories are ordered separately. For example Left and Right Cotsides are order individually and not as a pair. If unsure ask your Invacare Product Specialist or Customer Services.
- If purchasing an Accessory to add to existing beds, it is important to specify the model of Bed (e.g. Alegio) and if possible the Serial Number/Year the bed was purchased. Not all Accessories fit to all Beds. Customer Services will be able assist in identifying the correct Accessory for you. Not all Accessories fit all beds, the model and year can make a difference.
- Images on this brochure may differ from actual accessories supplied. This will depend on the bed it is to be fitted to. But the basic functionality and purpose of the accessories is as describe and is the same for all Beds.



### Local distributor



# Yes, you can.®

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Invacare New Zealand,  
PO Box 62-124, 4 Westfield Place, Mt Wellington, Auckland, New Zealand.  
Freephone: 0508 468 222, Freefax: 0508 807 788 Email: sales@invacare.co.nz, www.invacare.co.nz

All information is subject to change without prior notice. Photographs are for illustration purposes only and may vary from product description. Invacare® is a registered trademark and Yes, you can.™ is a trademark of Invacare® Corporation. 2005